noves BINO

Bedienungsanleitung
User Manual
Mode d'emploi
Istruzioni d'uso
Manual de instrucciones
Gebruiksaanwijzing
Betjeningsvejledning
Bruksanvisning
Bruksanvisning
Käyttöopas
Instrukcja obsługi
Provozní návod
Руководство по эксплуатации
サービスマニュアル

ESCHENBACH

Deutsch.....	3
English.....	4
Français.....	5
Italiano.....	6
Español.....	7
Nederlands.....	8
Dansk.....	9
Svenska.....	10
Norsk.....	11
Suomi.....	12
Polski.....	13
Česky.....	14
Русский.....	15
日本語.....	16

Deutsch

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English

Congratulations on your purchase of the novel BINO high magnification reading glasses from Eschenbach. The reading glasses novel BINO are supplied in a binocular design. The lenses are protected by a cera-tec® hard coating. Scratch-resistance is an important criterion for the quality and service life of the lens.

Important safety instructions, please read them attentively!

- ▶ **Risk of fire! Lenses in optical devices can cause significant damage through the "burning-glass effect" if handled or stored improperly. Take care to never to leave optical lenses lying in direct sunlight.**
- ▶ **Risks of glare and injury! Never look directly into the sun or other bright light sources when using optical devices.**
- ▶ **Also inform other people and children especially of this safety measure!**
- ▶ **Protect your reading glasses against impacts or shocks and excess heat. Never place your reading glasses on radiators or in direct sunlight.**
- ▶ **These reading glasses are suitable only for reading and may not be used in road traffic.**
- ▶ **The product is a Class I medical device (non-sterile, without measuring function) and complies with Regulation (EU) 2017/745. For use of the product outside the EU member states, the respective country-specific directives must be observed.**

Utilisation

1. Use appropriate lighting, one which allows you to evenly illuminate your reading material and does not cause reflections.
2. Put your reading glasses on.
3. To find the correct reading distance, hold the reading material close to your glasses.
4. Slowly move the reading material away from the glasses until the print is seen as a sharp image. This is the correct reading distance.
 - The print appears blurred when the reading distance is too large or too small.

Reading distances for different reading glasses

- The reading distance is dependent on the magnification of the reading glasses being used.

Article number	Diopters	Reading distance
16824 / 168204	4 dpt.	approx. 25 cm
16825 / 168205	5 dpt.	approx. 20 cm
16826 / 168206	6 dpt.	approx. 17 cm
16828 / 168208	8 dpt.	approx. 12,5 cm
168210/ 1682010	10 dpt.	approx. 10 cm

If you have questions, please always contact your dealer.

Care Instructions!

- ▶ **To clean the reading glasses, do not use alcohol or organic solvents.**
- ▶ **Do not clean the reading glasses under running water or in an ultrasonic bath.**

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Eschenbach Optik GmbH