
ridoMED/teleMED



Bedienungsanleitung
User Manual
Mode d'emploi
Istruzioni d'uso
Manual de instrucciones
Gebruiksaanwijzing
Betjeningsvejledning
Bruksanvisning
Käyttöopas
Instrukcja obsługi
Provozní návod
Руководство по эксплуатации
サービスマニュアル

 **ESCHENBACH**

Deutsch.....	4
English.....	6
Français.....	8
Italiano.....	10
Español.....	12
Nederlands.....	14
Dansk.....	16
Svenska.....	18
Norsk.....	20
Suomi.....	22
Polski.....	24
Česki.....	26
Русский.....	28
日本語.....	30

Deutsch

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English

You have selected a quality product from the house of Eschenbach which has been produced using the latest manufacturing processes in Germany, a top-class „Made in Germany“ product. We congratulate you on your choice.

Intended use

Various pre-assembled binocular systems are available from Eschenbach:

Magnifying spectacles for near viewing
ridoMED

— binocular (with convergence positions)

Magnifying spectacles for far viewing
teleMED

— binocular (with parallel positions)

Safety instructions

- ▶ **Danger of fire** Lenses used in optical equipment may cause considerable damage through their „burning glass effect“ if used or stored incorrectly. Ensure that optical lenses are never left uncovered in the sunlight!
- ▶ **Danger of blinding or injury** Never look at the sun through any optical equipment!
- ▶ **Danger of tripping when wearing magnifying spectacles or telescopic spectacles!**
- ▶ **Protect your magnifying spectacles or telescopic spectacles from knocks, bangs, moisture and excessive heat!**

- ▶ **When focusing the lenses, do not rotate past the limit stop!**
- ▶ **The product is a Class I medical device (non-sterile, without measuring function) and complies with Regulation (EU) 2017/745. For use of the product outside the EU member states, the respective country-specific directives must be observed.**

Adjusting the binocular system

Anatomical adjustments are carried out by the eyeglass optician.

Anatomical adjustments

1. Adjust side length.
2. Inclination (Frame angle): adjust to suit head position of user

Adjustment of pupil distance (PD)

1. Adjust the near/far PD using the scale [2] set on the eye-side of the binocular carrier [1].
 - The viewing field of both eyes must coincide, i.e. only **one** clear image must appear.
Fit the Galilei systems as closely as possible to the eye to ensure the maximum possible viewing field for the user.
2. Once the PD has been set, lock the eyepiece holders [3] with the clips [4].

Adjust for viewing errors

The system can be corrected at the eye-side (far correction).

- For eye-side correction, use the accessories from Article No. 16221 (Former for spectacle glasses with 4 replacement attachment rings).
- 1. Grind the correction lenses to 22 mm.
- 2. Attach correction lenses using the attachment rings to the eye-side of the system.
- If eye-side correction is not necessary, the eyepiece cups [5] can be fitted with caps [6].

If required, the adjustment supplement in the Eschenbach brochure "Adjustment of magnifying viewing aids" Article No. 162510 gives additional information.

- 3. The focusing ring can be rotated [7] for accurate focusing.
- The focusing positions in the system can be used to compensate spherical correction (Dioptric compensation ± 5 Dpt.).

The focusing ring is difficult to rotate in order to prevent the focus setting from changing by itself.

Instructions for care

The magnifying/telescopic spectacles should be kept in their case when not in use.

Clean the lenses with a soft, lint-free cloth (e. g. lens cleaning cloth), for harder to remove dirt (e. g. finger prints) slightly dampen the cloth.

Technical data

ridoMED	16362	16363	16364
Magnification	2.5 ×	3.0 ×	4.0 ×
Viewing field values *) Near (mm)	91/350	54/200	34/250
Viewing angle	15.0°	12.5°	8.0°
Working distance (mm)	≈ 350	≈ 200	≈ 250
Lens Ø	23 mm		
Weight	70 g		
Dioptric compensation	± 5 dpt.		

teleMED	1634	16344
Magnification	3 ×	4 ×
Viewing field values *) Far (mm)	160/1000	130/1000
Viewing angle	9.5°	7.5°
Working distance (mm)	700 ~ ∞	780 ~ ∞
Lens Ø	23 mm	
Weight	70 g	
Dioptric compensation	± 5 dpt.	

*) The effective viewing field in the Galilei systems is dependent on the distance between the eye and the eyepiece. Therefore, the viewing field values given are average values (Viewing angle was measured at HSA [cornea-lens distance] 10 mm).

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Warranty

Within the framework of the applicable statutory provisions, we provide a warranty for the functioning of the product described in this manual with regard to faults arising in connection with manufacturing errors or material errors. In event of damage due to improper handling, including damage from a fall or impact, no warranty claim exists. Warranty claims are only accepted upon presentation of the sale receipt!

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Eschenbach Optik GmbH